

A young man with short brown hair is captured in a dynamic running pose on the deck of a boat. He is wearing a black long-sleeved athletic top with grey panels on the chest and shoulders, black shorts, white socks with black stripes, and grey and white sneakers. The background shows a clear blue sky and the dark blue sea. The boat's railing is visible in the foreground.

POINT

**NAVIGATING
YOUR PATH TO
OPTIMAL**
Health



Nisantasi Orthopedics is a well-established company founded in 1993, specializing in the manufacturing and distribution of orthopedic products. We are among the leading healthcare manufacturers in Turkey, dedicated to providing high-quality and innovative solutions to our customers.

As **Nisantasi Orthopedics**, we own the brands “**Point**” and “**Orthez**”. We hold all necessary quality certifications, including ISO 13485, ISO 9001 and CE, ensuring compliance with international standards. Maintaining these standards is fundamental in delivering trust and quality to our customers.

In addition to serving as OEM manufacturers for numerous brands both domestically and internationally, we export our products to 20 countries worldwide. To support our international presence, we have established distributorships in many countries. We take pride in our strong global presence and strive to expand it further.

At **Nisantasi Orthopedics**, we prioritize customer satisfaction by understanding and meeting their needs and expectations. Our mission is to continue developing innovative and high-quality products to enhance the health and comfort of our customers. We aim to maintain a leading role in the healthcare sector by prioritizing customer satisfaction.

Nisantasi Orthopedics Products is committed to continuous improvement and remaining open to innovations to play a significant role in the healthcare sector.

CONTENTS

4-5



COLLAR

7-10



SHOULDER & ELBOW

12-17



WRIST

19-25



CORSET

27-30



KNEE & HIP

32-36




ANKLE

38-39



FOOT HEALTH

41-42



COMPRESSION STOCKING

44



OTHER PRODUCTS

46



PILLOW


48



WALKING AID



COLLAR



101 SOFT FOAM COLLAR

It is produced by cotton sheathed sponge fabric in accordance with the anatomical neck structure. With its elastic structure, it restricts neck movements, controls them and keeps the soft tissues loose by keeping the area warm.

INDICATIONS

- Used In soft tissue crushes
- Used In neck eclipses
- Used At the beginning of neck hernia
- Used for thermotherapy



S	M	L	XL
28-33	34-38	39-43	44-50
PERIMETER / CM			
COLOR <input checked="" type="radio"/> BLACK <input type="radio"/> BEIGE			



102 CERVICAL COLLAR WITH CHINS

It is made of plastic material, designed according to the anatomical shape of the neck. Limits neck movements with its chin support. Height adjustment can be made with the help of the Velcro straps on the top. It reduces pressure on the neck muscles by fixing the neck area.

INDICATIONS

- Used in neck hernia
- Used in neck calcifications
- It is used for post-operative protection and stabilization



S	M	L	XL	XXL
32-35	36-39	40-43	43-46	46-50
PERIMETER / CM				



103 CERVICAL COLLAR WITHOUT CHINS

It is made of plastic material, designed according to the anatomical shape of the neck. Height adjustment can be made with the help of the Velcro straps on the top. It reduces pressure on the neck muscles by fixing the neck area.

INDICATIONS

- Used in neck hernia
- Used in neck calcifications
- It is used for post-operative protection and stabilization



S	M	L	XL	XXL
32-35	36-39	40-43	43-46	46-50
PERIMETER / CM				



104 NELSON COLLAR

It is made of durable plastozote material with a thickness of 12 mm in accordance with the anatomical neck structure. It controls neck movements, especially by restricting it to the front and back, and keeps the soft tissues loose by keeping the area warm.

INDICATIONS

- Used in soft tissue crushes
- Used in neck eclipses
- Used in neck hernia
- Used for post-traumatic protection
- Used for thermotherapy



S	M	L	XL
28-33	34-38	39-43	44-50
PERIMETER / CM			

COLLAR



105 PHILADELPHIA COLLAR

It is produced from 12 mm thick plastozote material in accordance with the anatomical neck structure. The neck collar, which is applied to the neck in the form of two parts, completely restricts neck movements and keeps the soft tissues loose by keeping the area warm.

INDICATIONS

- Used in soft tissue crushes
- Used in vertebral fractures
- Used in neck calcifications
- Used in neck hernia
- First aid to trauma
- Used in post-surgical rehabilitation



S	M	L
33-37	38-42	42-46
PERIMETER / CM		



106 PHILEDELPHIA TRACHEOSTOMY COLLAR

It is produced from 12 mm thick plastozote material in accordance with the anatomical neck structure. The neck collar, which is applied to the neck in the form of two parts, completely restricts neck movements and keeps the soft tissues loose by keeping the area warm.

INDICATIONS

- Used in soft tissue crushes
- Used in vertebral fractures
- Used in neck calcifications
- Used in neck hernia
- First aid to trauma
- Used in post-surgical rehabilitation



S	M	L
33-37	38-42	42-46
PERIMETER / CM		



108 MINERVA ORTHOSIS

It is made of polyethylene plastic and cotton fabric material. The neck collar, which is applied to the neck and upper body in a multi-part manner, restricts the neck movements in the desired position by taking support from the shoulder and upper body elements. With the help of its adjustable mechanism, it helps apply the desired extension and flexion angle of the neck completely and safely

INDICATIONS

- Used in neck hernia
- Used in neck calcifications
- It is used for post-operative protection and stabilization



S	M	L
28-33	34-39	40-45
PERIMETER / CM		



111 FOAM COLLAR WITH PLASTIC SUPPORT

It is produced by cotton sheathed sponge fabric in accordance with the anatomical neck structure. With its elastic structure and extra plastic support, it restricts neck movements, controls them and keeps the soft tissues loose by keeping the area warm.

INDICATIONS

- Used in soft tissue crushes
- Used in neck eclipses
- Used at the beginning of neck hernia
- Used for thermotherapy



S	M	L	XL
28-33	34-38	39-43	44-50
PERIMETER / CM			

SHOULDER & ELBOW





152 ARM SLING (MESH)

It is made of breathable airtex mesh fabric and is designed in accordance with arm and shoulder anatomy. It enables the patient to keep his/her arm close to the body and parallel to the ground. Suitable for left and right arm use.

INDICATIONS

- Used in shoulder dislocations and semi-dislocations
- Treatment of forearm diseases
- Used for resting after surgery
- It is used to prevent edema formation in the fracture



CODE	S	M	L	XL	XXL
152	30-33	33-36	36-39	39-42	42-45
153	ONE SIZE				
PERIMETER / CM					



155 ARM SLING


It's inner surface is made of terry cloth and its outer surface is sponge fabric and it is designed in accordance with arm and shoulder anatomy. It ensures that the patient's arm is stable, parallel to the ground, and close to the body. Suitable for left and right arm use .

INDICATIONS

- Used in shoulder dislocations and semi-dislocations
- Treatment of forearm diseases
- Used for resting after surgery
- It is used to prevent edema formation in the fracture



S	M	L	XL	XXL
30-33	33-36	36-39	39-42	42-45
PERIMETER / CM				



156 ARM SLING (SIMPLE)

Its inner surface is made of terry cloth and its outer surface is sponge fabric and it is designed in accordance with arm and shoulder anatomy. With the help of its practical design, it allows the patient to keep his arm close to the body and parallel to the ground. Suitable for left and right arm use.

INDICATIONS

- Used in shoulder dislocations and semi-dislocations
- Treatment of forearm diseases
- Used for resting after surgery
- It is used to prevent edema formation in the fracture

**ONE
SIZE**




157 ARM SLING (30°)

It is a padded arm sling at a 30 degree angle. It is made of sponge fabric and With the help of its design, it allows the patient's arm to stay parallel to the ground close to the body. Suitable for left and right arm use. There is an additional belt for the waist circumference to facilitate daily use.

INDICATIONS

- Used in rotator cuff muscle tears
- Used in shoulder calcifications
- Used for rehabilitation after shoulder implant prosthesis
- Used for rehabilitation after frozen shoulder arthroscopy

**ONE
SIZE**



158 ARM SLING (45°-60°)

It is a padded arm sling at a 45-60 degree angle. It is made of sponge fabric and With the help of its design, it allows the patient's arm to stay parallel to the ground close to the body. Suitable for left and right arm use. There is an additional belt for the waist circumference to facilitate daily use.

INDICATIONS

- Used in rotator cuff muscle tears
- Used in shoulder calcifications
- Used for rehabilitation after shoulder implant prosthesis
- Used for rehabilitation after frozen shoulder arthroscopy

**ONE
SIZE**



159 CLAVICLE BANDAGE

It is made of sponge fabric with a Terry cloth on the inner surface. With the Velcro adjustment, it allows the shoulders to be pulled back with the desired tension and provides a smooth posture. With the help of its clavicle bandage design, it provides the necessary immobilization and prevents dislocation.

INDICATIONS

- Used in Clavicle (collarbone) fractures
- Used in dislocations of the joint between the clavicle (collarbone) and sternum (breastbone)
- Used in shoulder joint injuries and posture disorders (postural kyphosis)



CODE	XS	S	M	L	XL	XXL
159	18-24	25-30	31-36	37-42	43-48	49-54
160	ONE SIZE					
<i>PERIMETER / CM</i>						



161 VELPEAU BANDAGE

It is made of sponge fabric and is designed in accordance with the arm and shoulder anatomy. It provides immobilization of the shoulder and upper arm region in abduction. Bandage wrapping around the chest area enhances fixation. Velcro on the front provides ease of use.

INDICATIONS

- Used in Shoulder dislocations and semi-dislocations
- Used in tears of the rotator cuff muscles
- Used in Fractures of the humerus (arm bone) and scapula
- Used for the Soft tissue injuries in the shoulder and arm



CODE	S	M	L	XL	XXL
161	30-33	33-36	36-39	39-42	42-45
162	ONE SIZE				
<i>PERIMETER / CM</i>					



165 SHOULDER SUPPORT

It is made of neoprene material and designed in accordance with the anatomy of the shoulder. It is a shoulder bandage that protects the heat in the area and increases blood circulation thanks to the compression it applies to the shoulder area. Thanks to the applied pressure, it prevents edema. Suitable for use with both arms.

INDICATIONS

- Used in Shoulder bone fractions
- Used in Shoulder arthritics
- Used in post-operative irritations
- Used for rehabilitation after surgery

**ONE
SIZE**




167 HUMERUS BRACE

It is a shoulder orthosis made of sponge fabric with a Terry cloth on the inner surface. The plastic and aluminum supports inside fix and protect the humerus. These supports can be applied to the desired shoulder size. It is used to protect and stabilize the humerus in shoulder dislocations, after surgery and plaster cast. An optional arm hanger is available.

INDICATIONS

- Used for protection and stabilization in shoulder dislocations
- Used In humeral head fractures
- Used for rehabilitation after surgery
- Used for rehabilitation after shoulder plaster application



S	M	L
22-24	24-27	27-30
PERIMETER / CM		




168 CONTRACTOR ARM ORTHOSIS

It is made of thermoplastic material and is designed in accordance with shoulder and arm anatomy. With the help of the mechanism in the elbow area, the desired adjustment can be made. It can be adjusted in degrees of flexion (0-15-30-45-60-75-90-105) and extension (-15-0-15-30-45-60-75-90). With the aluminum support that prevents the wrist from falling, length adjustment can be made on the orthosis according to the arm length. Shoulder strap provides comfort in use.

INDICATIONS

- Used for Fractures of internal or external epicondyles
- Used After ligament injuries
- Used After elbow implant prosthesis application
- Used to speed up the healing process
- Used for conservative treatment after elbow dislocation



S	M	L
22-24	24-27	27-30
PERIMETER / CM		




169 EPICONDYLITIS ELBOW SUPPORT

It is made of neoprene material and is designed in accordance with the anatomy of the elbow. With help of the special 3-point silicone pad that provides local compression, it reduces the pain in the tendon. Designed for left and right arm use, it accelerates blood circulation and preserves body temperature.

INDICATIONS

- Used for Lateral epicondylitis
- Increased intra-articular fluid
- Used in painful conditions of the elbow



CODE	S	M	L	XL	XXL
169	22-24	24-26	26-28	28-30	30-32
170	ONE SIZE				
PERIMETER / CM					



174 ELBOW IMMOBILIZER

It is made of breathable neoprene fabric. It prevents sweating thanks to towel laminated inside surface. Thanks to the steel bars, it ensures that the arm stops in stabilization and the immobilization of the elbow in extension. It is anatomically designed and supportive after elbow injuries and surgery.

INDICATIONS

- Elbow joint traumas.
- Rest after surgery.
- It helps protection and immobilization.

**ONE
SIZE**



175 EPICONDYLITIS BANDAGE

It is made of neoprene fabric and designed in accordance with the anatomy of the elbow. Thanks to its special pad that provides local compression, it reduces the pain in the tendon. Designed for left and right arm use, it accelerates blood circulation and preserves body temperature. It is Universal size and can be adjusted individually.

INDICATIONS

- Tendinitis
- Golfer's elbow
- Tennis elbow
- It is used for post-operative protection and rest.

**ONE
SIZE**



226 KNITTED ELBOW SUPPORT

It is produced using flexible, breathable, 3D knitting technology and is designed in accordance with the anatomy of the elbow. It reduces tendon tension and pain with local compression.

INDICATIONS

- Increased intra-articular fluid
- Provides protection and support to the elbow
- Arthritis
- It is suitable for daily use against pain

S	M	L	XL	XXL
22-24	24-26	26-28	28-30	30-32
PERIMETER / CM				



227 KNITTED EPICONDYLITIS ELBOW SUPPORT


It is produced using flexible, breathable, 3D knitting technology and is designed in accordance with the anatomy of the elbow. It provides compression on the muscle groups attached to the lateral epicondyle

INDICATIONS

- In lateral epicondylitis
- When intra-articular fluid increases
- It is used in painful conditions of the elbow

S	M	L	XL	XXL
22-24	24-26	26-28	28-30	30-32
PERIMETER / CM				






200 WRIST SPLINT

It is made of breathable anti-allergic fabric and is designed in accordance with the anatomy of the wrist. It has an anatomically shaped aluminum support that provides immobilization of the wrist in neutral position. There is a bandage that increases the stabilization around the wrist.

INDICATIONS

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Used for resting in painful conditions of the wrist
- Used for injury and post-operative protection



S	M	L	XL	XXL
14-16	16-18	18-20	20-22	22-24
PERIMETER / CM				
COLOR <input checked="" type="radio"/> BLACK <input type="radio"/> BEIGE				




201 WRIST SPLINT WITH THUMB

It is made of breathable anti-allergic fabric. It is designed in accordance with the anatomy of the wrist and thumb. With its elastic structure, it strengthens the muscles by compressing the tissues and prevents swelling. With the help of aluminum splint, it completely restricts the movement of the wrist and thumb and takes it to rest. There is a bandage around the wrist that increases stabilization.

INDICATIONS

- Used When carpal tunnel and dequervain are together
- Used to Grade joint calcifications
- Used for Protection in thumb ligament laxity
- Used after plaster in fractures



S	M	L	XL	XXL
14-16	16-18	18-20	20-22	22-24
PERIMETER / CM				
COLOR <input checked="" type="radio"/> BLACK <input type="radio"/> BEIGE				



202 WRIST SPLINT MESH

It is made of breathable anti-allergic mesh fabric and is designed in accordance with the anatomy of the wrist. It has an anatomically shaped aluminum support that provides immobilization of the wrist in neutral position. There is a bandage around the wrist that increases stabilization.

INDICATIONS

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Used for resting in painful conditions of the wrist
- Used for injury and post-operative protection



S-M	L-XL
15-19	19-25
PERIMETER / CM	




203 WRIST SPLINT WITH THUMB MESH

It is made of breathable anti-allergic mesh fabric and is designed in accordance with the anatomy of the wrist. It has an anatomically shaped aluminum support that provides immobilization of the wrist in neutral position. There is a bandage around the wrist that increases stabilization.

INDICATIONS

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Used for resting in painful conditions of the wrist
- Used for injury and post-operative protection



S-M	L-XL
15-19	19-25
PERIMETER / CM	

HAND WRIST



204 WRIST SPLINT

It is made of breathable anti-allergic fabric and is designed in accordance with the anatomy of the wrist. It has an anatomically shaped aluminum support that provides immobilization of the wrist in neutral position. There is a bandage that increases the stabilization around the wrist.

INDICATIONS

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Used for resting in painful conditions of the wrist
- Used for injury and post-operative protection



**ONE
SIZE**



205 WRIST SPLINT WITH THUMB

It is made of breathable anti-allergic fabric and is designed in accordance with the anatomy of the wrist. It has an anatomically shaped aluminum support that provides immobilization of the wrist in neutral position. There is a bandage that increases the stabilization around the wrist.

INDICATIONS

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Used for resting in painful conditions of the wrist
- Used for injury and post-operative protection



**ONE
SIZE**



206 LACE WRIST SPLINT

It is made of airtex fabric with breathable mesh. It provides immobilization of the wrist in neutral position. There is an extra bandage on the wrist that increases stabilization. It is laced for ease of use.

INDICATIONS

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Used for resting in painful conditions of the wrist
- Used for injury and post-operative protection



**ONE
SIZE**



207 LACE WRIST SPLINT WITH THUMB


It is made of airtex fabric with breathable mesh. It provides immobilization of the wrist and thumb in neutral position. There is an extra bandage on the wrist that increases stabilization. It is laced for ease of use.

INDICATIONS

- Used When carpal tunnel and dequervain are together
- Used to Grade joint calcifications
- Used for Protection in thumb ligament laxity
- Used after plaster in fractures



**ONE
SIZE**



208 KNITTED WRIST SPLINT

It is produced using flexible, breathable, 3D knitting technology. The plastic support around the ulnar steloid reduces pressure on vessels and nerves while protecting the bony prominence. There is support on the volar face that provides immobilization of the wrist in neutral position. The bandage is designed for maximum fit, comfort and effective use, maintains body temperature, accelerates blood circulation and healing process.

INDICATIONS

- Carpal tunnel syndrome
- Tendinitis, tenosynovitis, tenovaginitis
- Painful conditions of the wrist
- Injury and post-surgery



S	M	L	XL	XXL
15-17	17-19	19-21	21-23	23-25
PERIMETER / CM				



210 DE QUERVAIN SPLINT

It is made of breathable anti-allergic fabric. It is designed in accordance with the anatomy of the hand, wrist and thumb. With its elastic structure, it strengthens the muscles by compressing the tissues and prevents swelling. With the help of aluminum splint, it fixes the thumb and restricts its movement.

INDICATIONS

- Used in Dequervain syndrome
- Used In grade joint calcifications
- Used for protection in thumb ligament laxity

**ONE
SIZE**



211 WRIST BANDAGE

It is made of breathable anti-allergic fabric. It is designed in accordance with the anatomy of the hand, wrist and thumb. With its elastic structure, it strengthens the muscles by compressing the tissues and prevents swelling. With the help of aluminum splint, it fixes the thumb and restricts its movement.

INDICATIONS

- Used In grade joint calcifications
- Used for protection in thumb ligament laxity

**ONE
SIZE**



211-1 WRIST BANDAGE WOVEN FABRIC

It is made of breathable anti-allergic woven fabric. It is designed in accordance with the anatomy of the hand, wrist and thumb. With its elastic structure, it strengthens the muscles by compressing the tissues and prevents swelling. With the help of aluminum splint, it fixes the thumb and restricts its movement.

INDICATIONS

- Used In grade joint calcifications
- Used for protection in thumb ligament laxity

**ONE
SIZE**

COLOR ● BLACK ● BEIGE



215 THE FROG SPLINT

Special aluminum finger splints with plastozote coated inner surface. It provides immobilization of the finger in extension. It is used for finger fractures, tendon and nerve injuries, strains and sprains.

INDICATIONS

- Used for Finger fractures
- Used for Tendon and nerve injuries
- Used for Finger injuries and sprains



S	M	L
4,5-5	5-6	6-7
PERIMETER / CM		



216 THE BASEBALL SPLINT

Special aluminum finger splints with plastozote coated inner surface. It provides immobilization of the finger in extension. It is used for finger fractures, tendon and nerve injuries, injuries and sprains.

INDICATIONS

- Used for Finger fractures
- Used for Tendon and nerve injuries
- Used for Finger injuries and sprains



S	M	L
4,5-5	5-6	6-7
PERIMETER / CM		



217 BUTONIER SPLINT

It is produced from thermoplastic material and designed in accordance with the finger anatomy. It is a finger orthosis that keeps the finger in extension according to the 3-point principle.

INDICATIONS

- Used for Correcting finger muscle immobility
- Used for Rehabilitation after surgery



S	M	L	XL
4	5	5-6	6-7
PERIMETER / CM			



218 MALLET FINGER SPLINT

It is made of special plastic material and is designed in accordance with finger anatomy. It is a finger orthosis that keeps the finger joint in extension. The finger is free from the PIP joint. It is used in extensor tendon injuries, fractures, sprains and sprains.

INDICATIONS

- Used in Mallet finger deformity
- Used in PIP joint fractures, strains, and sprains



1	2	3	4	5	6
4,5-5	5-5,5	5,5-6	6-6,5	6,5-7	7-8
PERIMETER / CM					




219 THERMOPLASTIC SPLINT

It is made of thermoplastic material and its inner surface is covered with terry cloth. Keeps the wrist and fingers in the resting position. It is used to stabilize the hand in the resting position after a plaster cast or in painful conditions of the wrist.

INDICATIONS

- Used After plaster cast in bone fractures joining the wrist
- Used For resting soft tissue injuries
- Carpal Tunnel Syndrome
- Used for Rehabilitation after surgery
- Used after Tendon and soft tissue damage



S	M	L
14-18	18-22	22-26
PERIMETER / CM		




221 THERMOPLASTIC ANTISPASTICITY SPLINT


It is made of thermoplastic material and its inner surface is covered with terry cloth. It is a hand and wrist orthosis made of thermoplastic material, the inner surface of which is covered with plastozote. Due to its corrugated finger design and resistance to tension, it is especially used in cases of spasticity in the hand.

INDICATIONS

- Spasticity
- Wrist discomforts



S	M	L
14-18	18-22	22-26
PERIMETER / CM		





222 THERMOPLASTIC DYNAMIC SPLINT


It is made of thermoplastic material and its inner surface is covered with terry cloth. With the help of its hard texture, it immobilizes the wrist in the anatomical position. With the help of aluminum and spring apparatus, open fingers do bending exercises.

INDICATIONS

- Used for Correcting muscle inactivity
- Used for correcting joint instability
- Used In ulnar nerve damage
- Used in post-surgical rehabilitation



S	M	L
14-18	18-22	22-26
PERIMETER / CM		





223 THERMOPLASTIC CLEINERT SPLINT


It is made of thermoplastic material and its inner surface is covered with terry cloth. The splint, designed in accordance with the anatomy of the wrist, allows passive flexion and active extension. It allows early controlled movement in flexor tendon injuries and after tendon repair.

INDICATIONS

- Used for Flexor tendon injuries
- Used for Rehabilitation after surgery



S	M	L
14-18	18-22	22-26
PERIMETER / CM		





224 MIDDLE FINGER FRACTURE SPLINT

Provides stabilization in the wrist and middle finger. It can be adjusted to the desired tightness with the velcro straps on it. There is an aluminum splint inside for maximum stabilization and immobilization.

INDICATIONS

- Used for Finger fractures
- Used for Tendon and nerve injuries
- Used for Finger injuries and sprains
- Used for Rehabilitation after surgery

**ONE
SIZE**



225 BOXER FRACTURE SPLINT

It provides stabilization in the wrist, little finger, and ring fingers. It can be adjusted to the desired tightness with the velcro straps on it. There are two aluminum splints inside for maximum stabilization and immobilization.

INDICATIONS

- MP joint fractures
- Used for Rehabilitation after surgery

**ONE
SIZE**

CORSET

POINT
NAVIGATING YOUR PATH TO OPTIMAL HEALTH



CORSET



250 CHEST CORSET

It is made of flexible, breathable, anti-perspiring, anti-allergic fishing line fabric. The chest corset supports the sternum area with terry cloth. There are two under wires on the chest. It supports and protects the sternum by reducing pain caused by sudden movement.

INDICATIONS

- Used In rib fractures
- Used for upper body stabilization
- Used after thoracic surgeries and traumas



CODE	S	M	L	XL	XXL
250	75-85	85-95	95-105	105-115	115-125
251	ONE SIZE				
PERIMETER / CM					



252 UMBILICAL HERNIA CORSET

It is made of flexible, breathable, anti-perspiring cotton fabric. It has a special designed pad that provides compression to the herniated area. It is applied in umbilical hernia.

INDICATIONS

- Used In abdominal soft tissue hernias
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation



CODE	S	M	L	XL	XXL
252	75-85	85-95	95-105	105-115	115-125
253	ONE SIZE				
PERIMETER / CM					



254 COLOSTOMY CORSET (26CM)

It is made of flexible, breathable, anti-perspiring, anti-allergic Elastic line fabric. It is used to reduce the risk of hernia after surgery in colostomy cases. It supports the colostomy bag and also the abdomen. It is produced with 26 cm.

INDICATIONS

- Protecting the colostomy bag
- Providing support after ostomy surgery
- Used to reduce the risk of parastomal hernia after surgery



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				



255 COLOSTOMY CORSET (16CM)


It is made of flexible, breathable, anti-perspiring, anti-allergic Elastic line fabric. It is used to reduce the risk of hernia after surgery in colostomy cases. It supports the colostomy bag and also the abdomen. It is produced with 16 cm.

INDICATIONS

- Protecting the colostomy bag
- Providing support after ostomy surgery
- Used to reduce the risk of parastomal hernia after surgery



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				



256 ABDOMINAL CORSET (16CM)

It is made of flexible, breathable, anti-perspiring, anti-allergic elastic line fabric. With its elastic structure, it protects and strengthens the muscles by compressing the tissues in the abdomen, and prevents swelling. It is produced with 16 cm.

INDICATIONS

- Used As a postpartum recovery
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation

CODE	S	M	L	XL	XXL
256	75-85	85-95	95-105	105-115	115-125
257	ONE SIZE				
PERIMETER / CM					




258 ABDOMINAL CORSET (26CM)

It is made of flexible, breathable, anti-perspiring, anti-allergic elastic line fabric. With its elastic structure, it protects and strengthens the muscles by compressing the tissues in the abdomen, and prevents swelling. It is produced with 26 cm.

INDICATIONS

- Used As a postpartum recovery
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation

CODE	S	M	L	XL	XXL
258	75-85	85-95	95-105	105-115	115-125
259	ONE SIZE				
PERIMETER / CM					



260 ELASTIC CORSET

With its elastic structure, it compresses the tissues, strengthens the muscles and maintains body temperature. It can be used underneath daily clothings.

INDICATIONS

- Used In abdominal organ cold
- Used in waist Calcification
- Used for low back pain
- Used for thermotherapy

S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				



261 ELASTIC CORSET WITH HOOK

With its elastic structure, it compresses the tissues, strengthens the muscles and maintains body temperature. It can be used underneath daily clothings.

INDICATIONS

- Used In abdominal organ cold
- Used in waist Calcification
- Used for low back pain
- Used for thermotherapy

S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				

CORSET



262 POSTUREX CORSET

It is made of flexible, breathable, non-sweating, anti-allergic elastic line fabric. It strengthens the muscles by providing compression to the tissues in the back and waist area with its elastic structure. With the help of plastic under wires on the back, it restricts the movement of the back vertebrae and puts it at rests.

INDICATIONS

- On the hunchback
- Used for Back muscle weakness
- Used in posture disorder



CODE	S	M	L	XL	XXL
262	75-85	85-95	95-105	105-115	115-125
263	ONE SIZE				
PERIMETER / CM					



264 MATERNITY CORSET

It is made of breathable and anti-perspiring elastic cotton material. Used before and after birth. This postpartum corset helps to recover large sagging in the anterior abdominal region. It gives a feeling of confidence and comfort for movement in daily life.

INDICATIONS

- Helps Back and waist pain during pregnancy
- Baby at risk of miscarriage
- Used in organ failure



CODE	S	M	L
264	70-110	90-130	115-145
265	ONE SIZE		
PERIMETER / CM			



266 LUMBOSACRAL CORSET WITH BELT (26CM)

It is made of flexible, breathable, anti-allergic, elastic 26 cm elastic line fabric. With its elastic structure, it strengthens the muscles by compressing the tissues in the waist and hips and prevents swelling. With the help of 4 steel under wires, it restricts the movement of the lumbar vertebrae and supports the vertebrae and puts the back at rest.

INDICATIONS

- Used In lumbar spine hernias
- Used in Sudden and long-term low back pain
- Used in post-surgical rehabilitation



CODE	S	M	L	XL	XXL
266	75-85	85-95	95-105	105-115	115-125
267	ONE SIZE				
PERIMETER / CM					



269 LUMBOSACRAL CORSET WITH BELT (32CM)


It is made of flexible, breathable, anti-allergic, elastic 32 cm elastic line fabric. With its elastic structure, it strengthens the muscles by compressing the tissues in the waist and hips and prevents swelling. With the help of 4 steel under wires, it restricts the movement of the lumbar vertebrae and supports the vertebrae and puts the back at rest.

INDICATIONS

- Used In lumbar spine hernias
- Used in Sudden and long-term low back pain
- Used in post-surgical rehabilitation



CODE	S	M	L	XL	XXL
269	75-85	85-95	95-105	105-115	115-125
270	ONE SIZE				
PERIMETER / CM					



271 DORSOLUMBAR CORSET

It is produced from flexible, breathable, anti-allergic elastic line fabric. With its elastic structure, it strengthens the muscles by compressing the tissues in the back, waist and hips. With the steel under wires on both sides of the spine, it helps to carry the upper body and rests by restricting the movement of the lumbar vertebrae appropriately.

INDICATIONS

- Used for spinal traumas
- Used in osteoporosis
- Used in Sudden and long-term back and waist pain
- Used in post-surgical rehabilitation



CODE	S	M	L	XL	XXL
271	75-85	85-95	95-105	105-115	115-125
272	ONE SIZE				
PERIMETER / CM					



273 NEOPRENE ABDOMINAL CORSET (26CM)

It is made of flexible, breathable neoprene fabric. With its elastic structure, it protects and strengthens the muscles by compressing the tissues in the abdomen, and prevents swelling. It is produced with 26 cm.

INDICATIONS

- Used As a postpartum recovery
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation

**ONE
SIZE**



274 T.L.S.O. CORSET

It is made of sponge fabric, the inner part of which is terry cloth, and there is a thermoplastic material that is portable and provides rigid support to the abdomen and hard part of the corset. The support in the sternum, which is different in this corset, which can be attached and removed like the Dorso-lumbar corset, provides flexion, extension, lateral tendency, full control and immobilization of the trunk in the post-op period. It can be extended or shortened according to the body length at the sternum.

INDICATIONS

- Used After severe trauma
- Used In vertebral fractures
- On the hunchback
- Used in post-surgical rehabilitation



S	M	L
60-85	85-110	110-140
PERIMETER / CM		



275 HYPEREXTENTION CORSET (JEWEEET)

It is a body orthosis that is made of special aluminum, the surfaces in contact with the skin are covered with terry cloth and according to the 3-point principle, it restricts the flexion of the body and encourages it in the hyperextension position. It limits the movements of the spine to a certain extent, with the help of its screw system on the orthosis, length and width adjustments can be made.

INDICATIONS

- Used for Fractures and traumas in the spine
- Used for kyphosis
- Used in post-operative rehabilitation



S	M	L
75-95	95-115	115-125
PERIMETER / CM		

CORSET



278 HERNIA TRUSS (DOUBLE)

It is designed according to the anatomical shape of the thigh region. It is made of tough elastic cotton fabric. There is a special support that applies pressure to the groin area. It is used in cases of unilateral inguinal hernia.

INDICATIONS

- Used in inguinal hernia



S	M	L
70-90	90-110	110-130
PERIMETER / CM		



279 HERNIA TRUSS (SINGLE)


It is designed according to the anatomical shape of the thigh region. It is made of tough elastic cotton fabric. There is a special support that applies pressure to the groin area. It is used in cases of unilateral inguinal hernia.

INDICATIONS

- Used in inguinal hernia



LEFT & RIGHT		
S	M	L
70-90	90-110	110-130
PERIMETER / CM		



280 POSTUREX CORSET COMFORT

It is made of flexible, breathable, non-sweating, anti-allergic elastic line fabric. It strengthens the muscles by providing compression to the tissues in the back and waist area with its elastic structure. With the help of plastic under wires on the back, it restricts the movement of the back vertebrae and puts it at rests.

INDICATIONS

- Used for Back muscle weakness
- Used in posture disorder

**ONE
SIZE**



282 LUMBOSACRAL CORSET WITH BELT COMFORT (26CM)


It is manufactured from breathable fabric. The height of corset is 26 cm. It supports waist and until to lower back. There are two plastic supports at the back of corset and two portable steel supports both sides of spine. 2 cross belts from the back to the navel area, provide additional stabilisation to the lumbar area. In this way, it relieves pain by reducing the load on the discs, and gives the patient a feeling of confidence and comfort.

INDICATIONS

- In herniated vertebral and calcification
- In sacralization and lumbalization
- Sacroiliac joint problems
- Sudden and prolonged low back pain
- Rehabilitation after surgical intervention



S-M	L-XL
75-100	100-125
PERIMETER / CM	



283 ABDOMINAL CORSET 3 BANDS

It is made of flexible, breathable, anti-perspiring, anti-allergic elastic line fabric. With its elastic structure, it protects and strengthens the muscles by compressing the tissues in the abdomen, and prevents swelling. It has 3 bants, all bants can apply the force in the region of use independently.

INDICATIONS

- Used As a postpartum recovery
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation



S-M	L-XL
75-100	100-120
PERIMETER / CM	
COLOR <input checked="" type="radio"/> BLACK <input type="radio"/> BEIGE	



284 ABDOMINAL CORSET 4 BANDS

It is made of flexible, breathable, anti-perspiring, anti-allergic elastic line fabric. With its elastic structure, it protects and strengthens the muscles by compressing the tissues in the abdomen, and prevents swelling. It has 4 bants, all bants can apply the force in the region of use independently.

INDICATIONS

- Used As a postpartum recovery
- Used for organ failure
- Supporting the abdominal wall
- Used in post-surgical rehabilitation



S-M	L-XL
75-100	100-120
PERIMETER / CM	
COLOR <input checked="" type="radio"/> BLACK <input type="radio"/> BEIGE	



285 KNITTED LUMBOSACRAL CORSET

It is produced using flexible, breathable, 3D knitting technology and is designed in accordance with the waist anatomy. There are 4 steel underwires at the back. Steel underwires can be tilted manually according to the lumbar lordosis. It adapts perfectly to the body with its conformity to the anatomical structure and superior quality. It reduces pain by reducing the tension in the muscles in the lumbar region.

INDICATIONS

- Lumbar spine hernias and calcifications
- In sacralization and lumbalization
- In sacroiliac joint problems
- Sudden and long-term low back pain
- In rehabilitation after surgery



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				



286 KNITTED LUMBOSACRAL CORSET WITH PAD

It is produced using flexible, breathable, 3D knitting technology and is designed in accordance with the waist anatomy. It is a corset with a specially designed triangular pad on the back lumbar part that creates a massage effect. Increases circulation with massage effect. It stabilizes the sacroiliac joint and the lumbosacral joint, and the compression applied reduces

INDICATIONS

- Lumbar spine hernias and calcifications
- In sacralization and lumbalization
- In sacroiliac joint problems
- Sudden and long-term low back pain
- In rehabilitation after surgery



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				

CORSET



290 ELASTIC TRICOT CORSET

With its elastic structure, it compresses the tissues, strengthens the muscles and maintains body temperature. It can be used underneath daily clothings.

INDICATIONS

- Used In abdominal organ cold
- Used in waist Calcification
- Used for low back pain
- Used for thermotherapy



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				



291 HIP CORSET

It is made of neoprene fabric. It is used as a support for warmth and compression in muscle tensions around the thigh and hip.

INDICATIONS

- Muscle strains around the thigh and hip
- Used as support for warmth and compression



S	M	L	XL	XXL
75-85	85-95	95-105	105-115	115-125
PERIMETER / CM				

COLOR BLACK BEIGE



281 WORKING CORSET (16CM)

It is 16 cm high corset manufactured from air fabric and supporting to waist and back. Consist 4 plastic supports at the back. 2 cross belts from the back to the navel area provides stabilization.

INDICATIONS

- In herniated vertebral and calcification
- In sacralization and lumbalization
- Sacroiliac joint problems
- Sudden and prolonged low back pain
- Rehabilitation after surgical intervention



S-M	L-XL-XXL
70-90	90-120
PERIMETER / CM	

KNEE & HIP

POINT
NAVIGATING YOUR PATH TO OPTIMAL HEALTH



KNEE & HIP



300 ELASTIC KNEE SUPPORT

It is made of acrylic fabric. With its elastic structure, it compresses the tissues, strengthens the muscles and maintains body temperature

INDICATIONS

- Used for knee calcification
- Used for knee pain
- Used for thermotherapy



S	M	L	XL	XXL
31-34	34-37	37-40	40-44	44-48
PERIMETER / CM				



301 KNEE SUPPORT CLOSED

It is made of neoprene fabric. It strengthens the muscles by compressing the tissues with its elastic structure and prevents swelling. It restricts the closing and opening movements of the knee and protects it from possible damage.

INDICATIONS

- Used for knee calcification
- Used In the increase of intra-articular fluid
- Used for Post-operative protection and rehabilitation
- Used For protection during sports activities
- Provides heat and compression



CODE	S	M	L	XL	XXL
301	31-34	34-37	37-40	40-44	44-48
302	ONE SIZE				
PERIMETER / CM					



305 KNEE SUPPORT WITH PATELLA

It is made of neoprene fabric. It is designed to support the knee joint. Patella is open around and supported by a silicone pad. With the help of the heat it spreads in the area and the compression it provides, it reduces the pain. It is suitable for use before and after surgery. It can be preferred for daily use as it does not put pressure on the kneecap.

INDICATIONS

- Used for knee calcification
- Used In the increase of intra-articular fluid
- Used for Knee cap pain
- Used for Post-operative protection and rehabilitation
- Used For protection during sports activities
- Provides heat and compression



CODE	S	M	L	XL	XXL
305	31-34	34-37	37-40	40-44	44-48
306	ONE SIZE				
PERIMETER / CM					



307 KNEE SUPPORT WITH PATELLA & LIGAMENT

It is made of neoprene fabric. It is designed to support the knee joint. Patella circumference is open and padded. Flexible under wires on the side increase stabilization and relieve the load on the knee. It further increases the stabilization with the Velcro tape on the thigh.

INDICATIONS

- Arthritis, bursitis, tendinitis
- Used In sprains and strains
- Used for meniscus problems
- Used for Ligament Instability
- For post-operative protection



CODE	S	M	L	XL	XXL
307	31-34	34-37	37-40	40-44	44-48
308	ONE SIZE				
PERIMETER / CM					



309 KNEE SUPPORT WITH PATELLA AND HINGES

It is made of neoprene material. Patella circumference is open and supported by a silicone pad. There are steel joints on both sides, providing maximum support and allowing controlled movement of the knee. It prevents knee hyperextension. It has Velcro straps that increase stabilization on the thigh and calf.

INDICATIONS

- Used for ligament injuries and ligament instability
- Used for meniscus problems
- Used In sprains and strains
- Restricts knee movement after surgery and speeds up the healing process



CODE	S	M	L	XL	XXL
309	31-34	34-37	37-40	40-44	44-48
310	ONE SIZE				
PERIMETER / CM					



311 KNEE SUPPORT WITH CROSSED LIGAMENTS

It is made of neoprene fabric. Designed to support knee joints, Patella circumference is open and supported with a silicone pad. With the help of semi-flexible supports on its inner and outer sides, it relieves the load on the lateral ligaments. It has an extra supporting bandage from the outside to the inside with the 3-point principle. It strengthens the muscles by compressing the tissues with its elastic structure, provides stabilization in the anterior cruciate ligaments and patella.

INDICATIONS

- Used in Anterior cruciate ligament injuries and knee instability
- Used for Patellar stabilization
- Used for meniscus problems
- Used In ligament damage
- Used for Post-operative protection



CODE	S	M	L	XL	XXL
311	31-34	34-37	37-40	40-44	44-48
326	ONE SIZE				
PERIMETER / CM					



314 PATELLAR BAND

It is made of neoprene material. It has an anatomically designed pad that supports the patellar tendon. It does not interfere with joint movements. The compression effect it applies reduces the tension and pain in the tendon. It can be used comfortably all day. It can be applied to both strings.

INDICATIONS

- Patellar tendinitis
- Os-good Schaltter disease
- Used for Patella femoral pain
- Chondromalacia patella

**ONE
SIZE**



315 KNITTED KNEE BRACE

It is produced using flexible, breathable, 3D knitting technology and is designed in accordance with the knee anatomy. The heat effect it provides on the knee joint provides long-term relief on the aching joints, muscles and surrounding tissues. Daily as it does not put pressure on the kneecap preferable for use.

INDICATIONS

- Arthritis, bursitis-increased intra-articular fluid
- Knee cap pain
- Meniscus problems
- First degree sprains and ligament instability
- Post-operative protection



S	M	L	XL	XXL
31-34	34-37	37-40	40-44	44-48
PERIMETER / CM				



316 KNITTED KNEE SUPPORT WITH PATELLA

It is produced using flexible, breathable, 3D knitting technology and is designed in accordance with the knee anatomy. The silicone pad around the patella relieves the pressure on the patella. Thanks to the heat it spreads in the area and the compression it provides, it reduces the pain. Ideal for people with kneecap pain when climbing, climbing or squatting stairs.

INDICATIONS

- Increased intra-articular fluid
- In rehabilitation after surgery
- Protection during and after sports activities
- Knee pain
- Arthritis



S	M	L	XL	XXL
31-34	34-37	37-40	40-44	44-48
PERIMETER / CM				



317 KNITTED KNEE SUPPORT WITH PATELLA & LIGAMENT

It is produced using flexible, breathable, 3D knitting technology and is designed in accordance with the knee anatomy. Patella is supported by a silicone pad. Flexible underwires on both sides increase stabilization and thus relieve the load on the knee. Thanks to the heat and compression it spreads in the area, it reduces the pain. It is suitable for use before and after surgery.

INDICATIONS

- Increased intra-articular fluid
- In rehabilitation after surgery
- In lateral ligament damage and kneecap pain
- Meniscus problems
- Arthritis, bursitis



S	M	L	XL	XXL
31-34	34-37	37-40	40-44	44-48
PERIMETER / CM				



318 KNEE IMMOBILIZER

It is made of sponge fabric with a terry cloth on the inner surface. It provides immobilization of the knee in extension. It has removable steel supports and Velcro that provide posterior, medial and lateral stabilization. Supportive for knee injuries and post-surgery, with the help of its three-piece pad design, it can be easily adjusted for leg circumference. It is available in 45 cm, 55 cm, and 65 cm lengths and is standard size.

INDICATIONS

- Used for Post-surgery in meniscus injuries
- Used for lateral ligament damage
- Used for Post-operative protection and immobilization
- Used for Rehabilitation after plaster cast



S	M	L
45	55	65
PERIMETER / CM		



319 KNEE BRACE WITH ADJUSTABLE ANGLE

It is an adjustable knee orthosis with an aluminum bar designed in accordance with the anatomical leg shape and terry cloth surfaces in contact with the skin. Flexion (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°) and extension (-15°, 0°, 15°, 30°, 45°, 60°, 75°, 90°). It provides excellent stabilization and support after knee injuries.


INDICATIONS

- Controlled post-operative injuries in ligament injuries and graded joint mobilization
- Post-surgery for meniscus injuries
- At rest at the desired angle after the surgical operation
- Used In knee recurvatum injuries



ADJUSTABLE

SHORT	LONG
47	58
PERIMETER / CM	



321 THIGH SUPPORT

It is made of neoprene fabric. Designed to support the thigh, with the help of the heat and compression it spreads to the area, it reduces the pain.

INDICATIONS

- Used in In muscle strains
- Used In painful conditions of the thigh
- For protection during and after sports activities
- Used to provide warmth and compression

**ONE
SIZE**



323 CALF SUPPORT

It is made of neoprene fabric. Designed to support the calf, with the help of the heat and compression it spreads to the area, it reduces the pain.

INDICATIONS

- Used in In muscle strains
- Used In painful conditions of the calf
- For protection during and after sports activities
- Used to provide warmth and compression

**ONE
SIZE**



370 HIP BRACE

It is a hip abduction orthosis that provides stabilization of the hip joint made of anatomically designed thermoplastic material. The inner surface in contact with the skin is covered with terry cloth. It can be adjusted up to 120 degrees of flexion in 15-degree steps. Abduction angle is adjustable.

INDICATIONS

- Used In hip calcification
- Used After total hip replacement or revision
- Used for stabilization purposes
- Used for Recurrent hip dislocation cases

**ONE
SIZE**



327 KNEE SUPPORT ELASTIC LONG

It is made of tricot fabric. With its elastic structure, it compresses the tissues, strengthens the muscles and maintains body temperature.

INDICATIONS

- Used for knee calcification
- Used for knee pain
- Used for thermotherapy




S	M	L	XL	XXL
31-34	34-37	37-40	40-44	44-48
PERIMETER / CM				



ANKLE





350 EIGHT BANDAGE

It is made of flexible, breathable, non-sweaty, anti-allergic, velcro-adhesive fabric, standard size and can be adjusted individually. Thanks to its elastic structure, it compresses the tissues in the foot, supports the muscles and relieves edema.

INDICATIONS

- Ankle sprains
- To support ligament weakness
- Mild soft tissue bruises
- Used for protection in sports activities

**ONE
SIZE**



351 ANKLE SUPPORT

It is made of neoprene fabric and is designed in accordance with the foot anatomy. It relieved the foot structure, relieved the edema and reduced the pain.

INDICATIONS

- Mild ligament injuries and instability
- Joint inflammations and calcifications
- For protective and support purposes after surgery and injury
- It is used in cases of edema, hematoma and effusion on the ankle



CODE	S	M	L	XL	XXL
351	21-23	23-25	25-27	27-29	29-31
352	ONE SIZE				
PERIMETER / CM					



353 ANKLE SUPPORT WITH MALLEOLAR PAD

It is made of neoprene fabric and is designed in accordance with the foot anatomy. There are pads on both sides of the ankle that cross the inner and outer malleolar. It relieved the foot structure, relieved the edema and reduced the pain.

INDICATIONS

- Mild ligament injuries and instability
- Joint inflammations and calcifications
- For protective and support purposes after surgery and injury
- It is used in cases of edema, hematoma and effusion on the ankle



CODE	S	M	L	XL	XXL
353	21-23	23-25	25-27	27-29	29-31
354	ONE SIZE				
PERIMETER / CM					



355 ANKLE SUPPORT WITH LIGAMENT

It is made of neoprene fabric and is designed in accordance with the foot anatomy. Anatomically designed plastic plates around the medial and lateral malleoli increase stabilization. It provides extra support with its elastic cross band, which is wrapped in the shape of eight. With the compression effect it applies, it supports the ankle and reduces edema and pain

INDICATIONS

- In ligament instability and tears
- In ankle sprains
- In post-plaster recovery
- Used for protection or support after surgery and injury



CODE	S	M	L	XL	XXL
355	21-23	23-25	25-27	27-29	29-31
356	ONE SIZE				
PERIMETER / CM					

ANKLE



357 ANKLE SUPPORT WITH PLASTIC STAY

Designed for ankle injury, rehabilitation and sports medicine. Plastic stay on both sides increases ankle stability. Reinforced crisscross strap pulls stay toward ankle to prevent sprain or swelling. Open heel applies less pressure on the Achilles tendon area.

INDICATIONS

- In ligament instability and tears
- In ankle sprains
- In post-plaster recovery
- Used for protection or support after surgery and injury



CODE	S	M	L	XL	XXL
357	21-23	23-25	25-27	27-29	29-31
358	ONE SIZE				
PERIMETER / CM					



359 KNITTED ANKLE SUPPORT

Flexible, breathable, produced using 3D knitting technology and designed in accordance with the foot anatomy. The heat effect it provides on the ankle provides long-term relief on aching joints, muscles and surrounding tissues. The compression effect it applies reduces edema and pain.

INDICATIONS

- In ankle sprains and strains
- To provide warmth and compression
- For compression purposes in post-injury edema



S	M	L	XL	XXL
21-23	23-25	25-27	27-29	29-31
PERIMETER / CM				



360 KNITTED ANKLE SUPPORT WITH MALLEOLAR PAD

Flexible, breathable, produced using 3D knitting technology and designed in accordance with the foot anatomy. There are silicone pads around the malleoli. By supporting the ankle, it reduces pain and edema with its compression effect. It is comfortable to use.

INDICATIONS

- In ankle sprains
- In rehabilitation after surgery
- In rehabilitation after plaster cast
- In ligament damage
- For compression purposes in post-injury edema



S	M	L	XL	XXL
21-23	23-25	25-27	27-29	29-31
PERIMETER / CM				



361 KNITTED ANKLE SUPPORT WITH ACHILLES TENDON

Flexible, breathable, produced using 3D knitting technology and designed in accordance with the foot anatomy. There is silicone support on both sides of the Achilles tendon and under the heel. It also has a portable heel placed under the other heel to prevent imbalance between both legs. It protects the heel and achilles tendon, reduces the tension in the tendon

INDICATIONS

- Achilles tendon pain
- After surgery or injury
- Retrocalcaneal bursitis
- Acute and chronic inflammations of the Achilles tendon



S	M	L	XL	XXL
21-23	23-25	25-27	27-29	29-31
PERIMETER / CM				



362 DORSI-FLEXION BANDAGE

It is made of neoprene fabric and is designed in accordance with the foot anatomy. It is an ankle brace used in low foot deformity that keeps the ankle in dorsiflexion at the desired tension. It provides a visible improvement in the patient's gait.

INDICATIONS

- Low foot deformity
- Peroneal neuropathy
- Muscle weakness in paralyzed patients
- Rehabilitation after surgery

ADJUSTABLE



363 AIRCAST ANKLE SUPPORT

The outer part is made of special thermoplastic material, and the inner part has pads with air cushions around the malleolus. It provides internal and external stabilization of the ankle.

INDICATIONS

- In recurrent ankle sprains in ligament damage
- Rehabilitation after surgery
- Protection during and after sports activities
- In rehabilitation after plaster cast

**ONE
SIZE**

COLOR BLACK WHITE



364 PADDED ANKLE SUPPORT

It is a supportive and comfortable ankle brace with a sponge pad designed in accordance with the joint shape. It provides internal and external stabilization in the ankle without hindering the natural joint movement. Ideal for the treatment and rehabilitation of mild to moderate ankle sprains.

INDICATIONS

- In recurrent ankle sprains in ligament damage
- Rehabilitation after surgery
- Protection during and after sports activities
- In rehabilitation after plaster cast

**ONE
SIZE**

COLOR BLACK WHITE



365 R.O.M WALKER

The base part is made of hard plastic material and the outer part is covered with a pile pad. It is an ankle orthosis that helps early mobilization with adjustable ankle dorsi and plantar flexion degrees. It is possible to adjust the plantar and dorsi flexion of the ankle at 0°-7.5°-15°-22.5°-30°-37.5°-45 degrees. The rocker model design on the sole supports the natural gait pattern, providing a more comfortable gait.

INDICATIONS

- It is used for postoperative stabilization in tendon, ligament and soft tissue injuries
- After acid tendon rupture surgery
- After heel and malleolar bone fractures



S	M	L
35-37	38-40	40-44
PERIMETER / CM		

ANKLE



366 ACHILLES BOOT

It is an ankle brace orthosis that is made of eva on the sole and hard plastic on the outside, covered with a terry pad inside, helping early mobilization. The plantar flexion height at the ankle is adjusted by the 3-layer pad on the heel.

INDICATIONS

- Chronic ligament instability
- Tendon, ligament and soft tissue injuries
- Achilles tendon tears for postoperative stabilization
- After surgery or injury



S	M	L
34-37	38-41	42-45
PERIMETER / CM		



367 REST AFO (FOOT GUARD)

It is made of thermoplastic material and is designed in accordance with the foot anatomy. It is a T-band ankle orthosis covered with terry cloth inside. By ensuring the ankle is fixed at 90 degrees, it prevents rotation with the T band.

INDICATIONS

- In chronic fasciitis syndrome
- Maintaining the foot position of the inpatient
- It is also used in rehabilitation after surgical intervention



S	M	L
35-37	38-40	40-44
PERIMETER / CM		



368 FOOT GUARD LUX (DEGREE ADJUSTABLE)

It is an ankle orthosis that provides gentle stretching of the ankle with dorsi flexion straps. The inner surface is covered with sponge fabric. Can be used for right and left feet. The non-slip surface on the sole provides ease of use for short distance walks.

INDICATIONS

- Acute and chronic inflammations of the Achilles tendon
- Plantar fasciitis
- Maintaining the position of the foot in the inpatient

ADJUSTABLE



S	M	L
35-37	38-40	40-44
PERIMETER / CM		



369 REFLEX AFO (RIGHT/LEFT)

It is an ankle orthosis made of plastic material and designed in accordance with the foot anatomy. It increases lateral stability by providing static dorsiflexion of the wrist. Velcro strap around the ankle provides comfort in use. It can be cut with scissors from the tip of the toe for a perfect fit into the shoe.

INDICATIONS

- Low foot deformity
- Muscle weakness in paralyzed patients
- Rehabilitation after surgery
- It is used in cases where it is desired to keep the ankle at 90 degrees



S	M	L
34-37	38-41	42-45
PERIMETER / CM		



421 HALLUX VALGUS SPLINT

It is a hallux valgus splint developed to prevent the progression of hallux valgus deformity. It consists of a thumb guard and elastic bandages covering the thumb. It prevents lateral deviation by applying a pulling force to the thumb. Elastic and soft material structure does not feel uncomfortable and provides maximum comfort. Suitable for shoes and daily use.

INDICATIONS

- Hallux valgus deformity

**ONE
SIZE**



422 HALLUX VALGUS NIGHT SPLINT

It is made of sponge fabric and velcro bands provide easy usage and stabilization. There is a support at the hallux valgus area provides correct position for toe. It is used in valgus deformity.

INDICATIONS

- Hallux valgus deformity
- Rehabilitation after surgical intervention

**ONE
SIZE**



423 HALLUX VALGUS SPLINT PLASTIC

It is a hallux-valgus support that adjusts flexion and extraction of the big toe. It can be easily attached and removed and can be used with shoes. The thumb is effectively supported. Healing after operation is accelerated.

INDICATIONS

AÇIKLAMA GELECEK

**ONE
SIZE**



FOOT HEALTH



400
SILICONE INSOLE



401
ARC HIGH SILICONE INSOLE



402
SILICONE INSOLE FABRIC COATED



403
MEDIAL AND LATERAL WEDGE



405
SILICONE HEEL CUP



406
SILICONE HEEL CUP FABRIC COATED WITH ADHESIVE



407
SILICONE TOE SEPARATOR



408
SILICONE ARCH SUPPORT WITH ADHESIVE



409
SILICONE BUNION PROTECTOR




FOOT HEALTH



411

SILICONE METATARSAL SUPPORT



414

SILICONE STRESS BALL



416

STARFLEX INSOLE



417

CORK INSOLE



418

CORK INSOLE WITH EPIN



419

STEEL INSOLE



424

HALLUX VALGUS SHOE SAĞ-SOL



425

PLASTER TOE



426

SABO SHOE

Compression Stockings

POINT
NAVIGATING YOUR PATH TO OPTIMAL HEALTH



DESCRIPTIONS

Varicose Socks manufactured from quality yarns and they are produced according to world standards and common sizes. It provides maximum comfort due to technological knitting technique and yarn quality and it provides maximum comfort and also it enables to minimize sweating to prevent bacterial growth thanks to respirable feature. Drop silicones, prevent slipping and does not disturb. Materials used in knitting are harmless to body and environment



551-555

CLASS I - CLASS II

VARICOSE STOCKING ABOVE KNEE

COLOR ● BLACK ● BEIGE



552-554

CLASS I - CLASS II

VARICOSE STOCKING PANTYHOSE

COLOR ● BLACK ● BEIGE



556

CLASS I - CLASS II

VARICOSE STOCKING PREGNANCY

COLOR ● BLACK ● BEIGE


550-553

CLASS I - CLASS II

VARICOSE STOCKING BELOW KNEE

COLOR ● BLACK ● BEIGE

Varicose
STOCKINGS



Leg Circumference	SIZES						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cF	36-46	37-51	42-57	45-60	49-65	52-68	55-70
cG	43-58	46-61	49-64	52-68	46-71	60-74	64-79

Antiembolism STOCKINGS




557
ANTI-EMBOLISM STOCKING
BELOW KNEE



558
ANTI EMBOLISM STOCKING
ABOVE KNEE



559
ANTI EMBOLISM STOCKING
BELTED



Leg Circumference	SIZES						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cF	36-46	37-51	42-57	45-60	49-65	52-68	55-70
cG	43-58	46-61	49-64	52-68	46-71	60-74	64-79



OTHER PRODUCTS

OTHER PRODUCTS



450
SILICONE BREAST PROTHESIS



281
NOSE GUARD MASK



371
SUSPENSORY SLIP



373
FREJKA PILLOW



375
PAVLIC HARNESS



452
TOILET SEAT RISER



700
CHAMBER INHALER



374
DKC ORTHOSIS




376
DENNIS BROWN SPLINT





PILLOW

PILLOW



600

VISCO NECK PILLOW



601

VISCO NECK GEL PILLOW



602

VISCO TRAVEL PILLOW



603

VISCO RING SEAT CUSHION
607-FOAM RING SEAT CUSHION




604

VISCO BACK CUSHION
609-FOAM VISCO BACK CUSHION



605

VISCO ANATOMIC SEAT CUSHION



606

REFLUX PILLOW



608

RING SEAT CUSHION (OPEN)



610

POSITION PILLOW



WALKING AID




650
WALKING STICK



651
TRIPOD WALKING STICK



653
CRUTCHES



654
CRUTCHES LUX



655
CANEDIAN CRUTCHES



656
WALKER (IRON)



657
WALKER (ALUMINUM)



658
FOLDING COMODE



661
MANUAL WHEELCHAIR

POINT

 **NIŞANTAŞI**
Ortopedi

NIŞANTAŞI ORTOPEDİ İML. İTH. İHR. SAN. VE TİC. LTD. ŞTİ.

19 Mayıs Mah. Saadet Hanım Sok. No: 13/A Şişli / İstanbul / Türkiye
Tel: +90 (212) 232 28 79 - +90 (212) 356 11 74 - +90 553 227 06 27
www.pointortopedi.com - info@pointortopedi.com